

# Anti-FOXD3 (5C4) Mouse antibody



Catalog Number: 168274

## Product name

Mouse IgG1

Anti-FOXD3 (5C4) Mouse antibody

## Applications

## Specificity

WB

Human, Mouse, Monkey

## Dilutions

## Antibody description

WB: 1/500 - 1/2000

FOXD3 (5C4) Mouse monoclonal antibody

ELISA: 1/10000

## Preparation

Antigen: Purified recombinant fragment of human FOXD3 expressed in E. Coli.

## Formulation

Ascitic fluid containing 0.03% sodium azide.

## Storage


Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

## Validations



Western blot detection of FoxD3 in Jurkat and HeLa cell lysates using FoxD3 mouse mAb (1:1000 diluted). Predicted band size: 48 kDa. Observed band size: 48 kDa.