

Product name

ASUN antibody

IP: 1:1000-1:10000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ASUN Rabbit Polyclonal antibody. Positive IP detected in MCF-7 cells. Positive WB detected in MCF7 cells. Positive IHC detected in human cervical cancer tissue, human testis tissue. Observed molecular weight by Western-blot: 70-80 kDa

Preparation

This antibody was obtained by immunization of ASUN recombinant protein (Accession Number: XM_017019632). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC, IP, WB

Dilutions


Recommended Dilution:

WB: 1:1000-1:10000


Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:108274(C12orf11 antibody) at dilution of 1:4000



Immunohistochemical of paraffin-embedded human cervical cancer using Catalog No:108274(C12orf11 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded

Catalog Number: 108274

human cervical cancer using Catalog No:108274(C12orf11 antibody) at dilution of 1:100 (under 40x lens)

IP Result of anti-C12orf11 (IP:Catalog No:108274, 4ug; Detection:Catalog No:108274 1:4000) with MCF-7 cells lysate 1040ug.

