

CHPF-N terminal antibody



Catalog Number: 109274

Product name

CHPF-N terminal antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CHPF-N terminal Rabbit Polyclonal antibody.
Positive WB detected in human brain tissue.
Observed molecular weight by Western-blot: 85-90 kDa, 65-70 kDa

Preparation

This antibody was obtained by immunization of CHPF-N terminal recombinant protein (Accession Number: BC008878). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109274(CHPF antibody) at dilution of 1:800