

MR1 Polyclonal Antibody



Catalog Number: 165674

Product name

MR1 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to MR1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 20-260 of human MR1 (NP_001181928.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB

Dilutions

WB 1:500 - 1:2000

Validations

Mouse skeletal muscle



Western blot - MR1 Polyclonal Antibody

Western blot analysis of extracts of mouse skeletal muscle, using MR1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.