

# C1orf226 antibody

Catalog Number: 108674

## Product name

C1orf226 antibody

## Specificity

Human; other species not tested.

## Antibody description

C1orf226 Rabbit Polyclonal antibody. Positive IP detected in A431 cells. Observed molecular weight by Western-blot: 28 kDa

## Preparation

This antibody was obtained by immunization of C1orf226 recombinant protein (Accession Number: NM\_001085375). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IP

## Dilutions

Recommended Dilution:

IP: 1:200-1:1000

## Validations



IP Result of anti-C1orf226 (IP:Catalog No:108674, 3ug; Detection:Catalog No:108674 1:300) with A431 cells lysate 1520ug.