

H9N2 Neuraminidase (His Tag) recombinant protein



Catalog Number: 501374

General Information

Protein Construction

A DNA sequence encoding the Influenza A virus (A/Hong Kong/1073/99(H9N2)) neuraminidase (NP_859038.1) (His36-Ile469) was expressed with a N-terminal polyhistidine tag.

Organism

H9N2

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant neuraminidase of Influenza A virus (A/Hong Kong/1073/99(H9N2)) comprises 451 amino acids and has a predicted molecular

mass of 50 kDa. The apparent molecular mass of the protein is approximately 65.9 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H9N2 (A/Hong Kong/1073/99) Neuraminidase / NA (His Tag) SDS-PAGE