

TSTA3 antibody

Catalog Number: 116374

Product name

TSTA3 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

TSTA3 Rabbit Polyclonal antibody. Positive WB detected in human colon tissue, HEK-293 cells, HeLa cells, human brain tissue, human kidney tissue. Observed molecular weight by Western-blot: 40kd, 68kd

Preparation

This antibody was obtained by immunization of TSTA3 recombinant protein (Accession Number: XM_011517269). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:116374(TSTA3 antibody) at dilution of 1:500