

BBS7 antibody

Catalog Number: 108374

Product name

BBS7 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

BBS7 Rabbit Polyclonal antibody. Positive IF detected in hTERT-RPE1 cells.

Preparation

This antibody was obtained by immunization of BBS7 recombinant protein (Accession Number: NM_018190). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IF

Dilutions


Recommended Dilution:

IF: 1:20-1:200

Validations



IF result (pericentriolar staining) of anti-BBS7 (Catalog No:108374, 1:50) with serum-starved hTERT-RPE1 by Dr. Moshe Kim.



IF result (along the length of cilia stain) of anti-BBS7 (Catalog No:108374) with serum-starved hTERT-RPE1 cells (1% PFA fixed) by Dr. Moshe Kim.