

ROBO3 antibody

Catalog Number: 114774

Product name

ROBO3 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ROBO3 Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue. Positive WB detected in mouse stomach tissue. Observed molecular weight by Western-blot: 70kd

Preparation

This antibody was obtained by immunization of ROBO3 recombinant protein (Accession Number: BC008623). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse stomach tissue were subjected to SDS PAGE followed by western blot with Catalog No:114774(ROBO3 antibody) at dilution of 1:200



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:114774(ROBO3 antibody) at dilution of 1:50 (under 10x lens)