

# Anti-FASLG antibody



Catalog Number: 141084

## Product name

Anti-FASLG antibody

## Specificity

Human Fas Ligand / FASLG / CD95L

## Antibody description

Rabbit Polyclonal to Human FASLG

## Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the C-terminus of the Human Fas Ligand / FASLG / CD95L, and purified by antigen affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, ICC/IF, IF

## Dilutions

WB: 5-20 µg/ml

ICC/IF: 0.5-3µg/mL

## Validations



Human Fas Ligand / FASLG / CD95L Immunofluorescence(IF) 15730



Lanes	A
Items	
Sample (whole cell lysate)	HeLa
Sample Volume (µg/lane)	30
Gel	13% SDS-PAGE reducing gel
Recommended Concentration	5-20 µg/ml
Secondary Antibody	Dylight 800-labeled Antibody To Rabbit IgG (H+L), at 1:5000 dilution.
Developed using Odyssey imaging system.	
Explanation	Predicted band size : 31 kDa Observed band size : 35 kDa

Human Fas Ligand / FASLG / CD95L Western blot (WB) 15729