

Anti-SFN antibody



Catalog Number: 142084

Product name

Anti-SFN antibody

Immunogen

[Human 14-3-3 sigma/Stratifin/YWHAS recombinant protein](#)

Specificity

Human 14-3-3 sigma / Stratifin / YWHAS

Antibody description

Mouse Monoclonal to Human SFN

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human 14-3-3 sigma / Stratifin / YWHAS . The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

FCM

Dilutions

FCM: 0.5-2 µg/Test

Validations



14-3-3 sigma / Stratifin / YWHAS Antibody, Mouse MAb, Flow Cytometry