

PPARD antibody

Catalog Number: 114084

Product name

PPARD antibody

ELISA, WB, IHC

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

IHC: 1:10-1:100


Specificity

human, mouse, rat; other species not tested.

Antibody description

PPARD Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue, COLO 320 cells, mouse skeletal muscle tissue. Positive IHC detected in human ovary tumor. Observed molecular weight by Western-blot: 54kd

Validations




Preparation

This antibody was obtained by immunization of PPARD recombinant protein (Accession Number: NM_177435). Purification method: Antigen affinity purified.

mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:114084(PPARD antibody) at dilution of 1:500

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:114084(PPARD antibody) at dilution of 1:10 (under 10x lens)

Applications