

Anti-CaMKII alpha (5D9) Mouse antibody



Catalog Number: 168084

Product name

Anti-CaMKII alpha (5D9) Mouse antibody

Specificity

Human, Mouse

Antibody description

CaMKII alpha (5D9) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human CaMKII- α protein fragments expressed in E.coli.

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of CaMKII- α in Mouse Brain lysates using CaMKII- α mouse mAb(dilution 1:1000).Predicted band size:50kDa.Observed band size:50kDa.