

HMGCR antibody

Catalog Number: 111484

Product name

HMGCR antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HMGCR Rabbit Polyclonal antibody. Positive IHC detected in human liver cancer tissue, human heart tissue.

Preparation

This antibody was obtained by immunization of HMGCR recombinant protein (Accession Number: NM_001130996). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:111484(HMGCR Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:111484(HMGCR Antibody) at dilution of 1:50 (under 40x lens)