

RASSF7 antibody



Catalog Number: 114484

Product name

RASSF7 antibody

Immunogen

[Human RASSF7 Recombinant protein \(GST tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

RASSF7 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, HeLa cells. Observed molecular weight by Western-blot: 36kd,40kd

Preparation

This antibody was obtained by immunization of RASSF7 recombinant protein (Accession Number: XM_017018361). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:114484(RASSF7 antibody) at dilution of 1:200