

# ARID3A Polyclonal Antibody



Catalog Number: 165484

## Product name

ARID3A Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to ARID3A

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 324-593 of human ARID3A (NP\_005215.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IP


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IP 1:20 - 1:50

## Validations




Western blot - ARID3A Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ARID3A antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - ARID3A Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat ileum using ARID3A antibody at dilution of 1:100



(40x lens).


# ARID3A Polyclonal Antibody



Catalog Number: 165484


Immunohistochemistry - ARID3A Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human oophoroma using ARID3A antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - ARID3A Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung cancer using ARID3A antibody at



dilution of 1:100 (40x lens).

Immunoprecipitation - ARID3A Polyclonal Antibody

Immunoprecipitation analysis of 150ug extracts of MCF7 cells using 3ug ARID3A antibody . Western blot was performed from the immunoprecipitate using ARID3A antibody at a dilution of 1:500.