

# Pre-IL18 antibody



Catalog Number: 107484

## Product name

Pre-IL18 antibody

## Specificity

Human, hamster; other species not tested.

## Antibody description

Pre-IL18 Mouse Monoclonal antibody. Positive WB detected in HeLa cells, HepG2 cells. Positive IF detected in MCF-7 cells, Hela cells. Observed molecular weight by Western-blot: 24 kDa

## Preparation

This antibody was obtained by immunization of Pre-IL18 recombinant protein (Accession Number: NM\_001386420). Purification method: Caprylic acid/ammonium sulfate precipitation.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications

ELISA, IF, WB


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100


## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107484(Pre-IL18 Antibody) at dilution of 1:1000



Immunofluorescent analysis of MCF-7 cells, using IL18 antibody Catalog No: at 1:25 dilution and Rhodamine-labeled goat anti-mouse IgG (red).



Immunofluorescent analysis of HeLa cells, using IL18 antibody Catalog No: at 1:25 dilution and Rhodamine-labeled goat anti-mouse IgG (red).