

# Anti-KAT5/Tip60 Mouse antibody



Catalog Number: 168484

## Product name

Mouse IgG1

Anti-KAT5/Tip60 Mouse antibody

## Applications

## Specificity

WB, IHC-P, ICC/IF, FC

Human

## Dilutions

## Antibody description

WB: 1/500 - 1/2000

KAT5/Tip60 Mouse monoclonal antibody

FC: 1/200 - 1/400

## Preparation


ELISA: 1/10000

Antigen: Purified recombinant fragment of human TIP60 (AA: 18-208) expressed in E. Coli.

## Validations

## Formulation

Purified antibody in PBS with 0.05% sodium azide



## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of Tip60 in MCF7,C6 and HeLa cell lysates using Tip60 mouse mAb (1:2000 diluted).Predicted band size:60KDa.Observed band size:60KDa.

## Clonality

Monoclonal

## Ig Type