

# Human DYNLL1/DLC1/DNCL1/DNCLC1/HDLC1 (His tag) recombinant protein



Catalog Number: 518484

## General Information

### Gene Name Synonym

Dynein Light Chain 1 Cytoplasmic, 8 kDa Dynein Light Chain, DLC8, Dynein Light Chain LC8-Type 1, Protein Inhibitor of Neuronal Nitric Oxide Synthase, PIN, DYNLL1, DLC1, DNCL1, DNCLC1, HDLC1

### Protein Construction

A DNA sequence encoding the Human DYNLL1/DLC1/DNCL1/DNCLC1/HDLC1 Met1-Gly89 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

### Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted

samples are stable at < -20°C for 3 months.

### Formulation

Lyophilized from a 0.2 µm filtered solution of 20mM TrisHCl, 200mM NaCl, 1mM DTT, pH 8.0.

## Usage Guide


### Storage

Lyophilized protein should be stored at <-20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days.

### Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in ddH<sub>2</sub>O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

## SDS-PAGE



Human DYNLL1/DLC1/DNCL1/DNCLC1/HDLC1 (His tag) recombinant protein