

Anti-Tryptophanyl tRNA synthetase/WRS Rabbit antibody



Catalog Number: 167884

Product name

Rabbit IgG

Anti-Tryptophanyl tRNA synthetase/WRS Rabbit antibody

Applications

WB, IHC-P, FC

Specificity

Human, Mouse, Rat

Dilutions

WB: 1/1000

Antibody description

Tryptophanyl tRNA synthetase/WRS Rabbit monoclonal antibody

IHC: 1/100-1/500

FC: 1/20


Preparation

Antigen: A synthetic peptide of human Tryptophanyl tRNA synthetase/WRS

Validations

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA



Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Western blot detection of Tryptophanyl tRNA synthetase/WRS in HeLa,A549,HL-60 cell lysates using Tryptophanyl tRNA synthetase/WRS Rabbit mAb(1:1000 diluted).Predicted band size:53kDa.Observed band size:53kDa.

Ig Type