

# EIF2B4 antibody

Catalog Number: 110184

## Product name

EIF2B4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

EIF2B4 Rabbit Polyclonal antibody. Positive WB detected in human skeletal muscle tissue, HeLa cells, HL-60 cells, K-562 cells, MCF7 cells.

Positive IHC detected in human breast cancer tissue, human testis tissue. Observed molecular weight by Western-blot: 60kd

## Preparation

This antibody was obtained by immunization of EIF2B4 recombinant protein (Accession Number: NM\_001034116). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:110184(EIF2B4 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:110184(EIF2B4 antibody) at dilution of 1:100 (under 10x lens)