

SH2D5 antibody

Catalog Number: 115184

Product name

SH2D5 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

SH2D5 Rabbit Polyclonal antibody. Positive IF detected in HEK-293 cells. Positive WB detected in HEK-293 cells, mouse kidney tissue. Observed molecular weight by Western-blot: 43-47 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001103160). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:20-1:200

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:115184(SH2D5 Antibody) at dilution of 1:300



Immunofluorescent analysis of HEK-293 cells using Catalog No:115184(SH2D5 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)