

# CXCL16 antibody



Catalog Number: 107184

## Product name

CXCL16 antibody

## Immunogen

[Human CXCL16 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

CXCL16 Mouse Monoclonal antibody. Positive WB detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 35 kDa, 60 kDa

## Preparation

This antibody was obtained by immunization of CXCL16 recombinant protein (Accession Number: NM\_022059). Purification method: Caprylic acid/ammonium sulfate precipitation.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:1000

## Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:107184(CXCL16 antibody) at dilution of 1:300