

ARHGEF16 antibody



Catalog Number: 108184

Product name

ARHGEF16 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ARHGEF16 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue. Positive WB detected in mouse kidney tissue. Observed molecular weight by Western-blot: 80 kDa

Preparation

This antibody was obtained by immunization of ARHGEF16 recombinant protein (Accession Number: XM_024446454). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:108184(ARHGEF16 antibody) at dilution of



1:100

Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:108184(ARHGEF16 antibody) at dilution of 1:50 (under 25x lens)