

Mouse sFRP4 (His Tag) recombinant protein



Catalog Number: 500584

General Information

Protein Construction

A DNA sequence encoding the mouse sFRP4 (NP_057896.1) (Met 1-Ser 351) was expressed with a C-terminal polyhistidine tag.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 97 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ala 22

Molecular Mass

The recombinant mouse sFRP4 consists of 341 amino acids after removal of the signal peptide and has a predicted molecular mass of 39.4 kDa.

In SDS-PAGE under reducing conditions, it migrates with an apparent molecular mass of 55-60 kDa due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse sFRP4 Protein (His Tag) SDS-PAGE