

USMG5 antibody



Catalog Number: 116584

Product name

USMG5 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

USMG5 Rabbit Polyclonal antibody. Positive IHC detected in human placenta. Positive WB detected in human liver tissue, human brain tissue. Observed molecular weight by Western-blot: 6kd

Preparation

This antibody was obtained by immunization of USMG5 recombinant protein (Accession Number: NM_001206427). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:116584(USMG5 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human placenta using Catalog No:116584(USMG5 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human placenta using Catalog No:116584(USMG5 antibody) at dilution of 1:100 (under 40x lens)