

Anti-Microsomal GST1 Rabbit antibody



Catalog Number: 168584

Product name

Rabbit IgG

Anti-Microsomal GST1 Rabbit antibody

Applications

Specificity

WB, IHC-P, FC

Human, Mouse, Rat

Dilutions

Antibody description

1:1000-1:5000

Microsomal GST1 Rabbit monoclonal antibody

1:20-1:50

Preparation

1:20

Antigen: A synthetic peptide of human Microsomal Glutathione S-transferase 1

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of Microsomal GST1 in Rat Liver cell lysates using Microsomal GST1 Rabbit mAb(1:1000 diluted).Predicted band size:18kDa.Observed band size:18kDa.