

MYO1D antibody

Catalog Number: 112984

Product name

MYO1D antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

MYO1D Rabbit Polyclonal antibody. Positive IP detected in mouse lung tissue. Positive WB detected in mouse lung tissue, mouse brain tissue. Observed molecular weight by Western-blot: 100 kDa

Preparation

This antibody was obtained by immunization of MYO1D recombinant protein (Accession Number: BC030602). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:112984(MYO1D antibody) at dilution of 1:200



IP Result of anti-MYO1D (IP:Catalog No:112984, 4ug; Detection:Catalog No:112984 1:300) with mouse lung tissue lysate 4000ug.