

XDH antibody

Catalog Number: 116984

Product name

XDH antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

XDH Rabbit Polyclonal antibody. Positive IHC detected in human heart tissue. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 147-150kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_000379). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:

WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:116984(XDH antibody) at dilution of 1:800