

SIAH1 antibody



Catalog Number: 115284

Product name

SIAH1 antibody

Specificity

human, mouse; other species not tested.

Antibody description

SIAH1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, A549 cells, COLO 320 cells, HepG2 cells, Jurkat cells, mouse small intestine tissue, RAW 264.7 cells. Observed molecular weight by Western-blot: 60 kDa

Preparation

This antibody was obtained by immunization of SIAH1 recombinant protein (Accession Number: XM_024450395). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115284(SIAH1 antibody) at dilution of 1:300