

# Anti-CNDP2 antibody



Catalog Number: 106284

## Product name

Anti-CNDP2 antibody

## Immunogen

[Mouse CNDP2 \(His Tag\) recombinant protein](#)

## Specificity

Mouse CNDP2

**Has cross-reactivity** in ELISA with Human CNDP2

## Antibody description

Rabbit monoclonal to CNDP2

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse CNDP2 (rM CNDP2; NP\_075638.2; Met1-Asn475).

## Formulation

0.2 µm filtered solution in PBS with 5% trehalose

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions

IHC-P: 1-10 µg/mL

## Validations



CNDP2 Antibody, Rabbit MAb,  
Immunohistochemistry

Immunochemical staining of mouse CNDP2 in mouse kidney with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections). Positive staining was localized to nephric tubule.