

Anti-RBAK antibody

Catalog Number: 176284

Product name

Anti-RBAK antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to RBAK

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human rbak. aa range:1-50

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations



Fig1.; Western Blot analysis of various cells using RBAK Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA).