

MLN51/CASC3 antibody



Catalog Number: 112684

Product name

MLN51/CASC3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MLN51/CASC3 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in A2780 cells, Jurkat cells, MCF7 cells, SKOV-3 cells. Observed molecular weight by Western-blot: 115-120 kDa

Preparation

This antibody was obtained by immunization of MLN51/CASC3 recombinant protein (Accession Number: XM_005257163). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



A2780 cells were subjected to SDS PAGE followed by western blot with Catalog No:112684(CASC3 antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells, using CASC3 antibody Catalog No:112684 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).