

PCBD2 antibody



Catalog Number: 113684

Product name

PCBD2 antibody

Immunogen

[Human PCBD2 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

PCBD2 Rabbit Polyclonal antibody. Positive WB detected in human skeletal muscle tissue. Observed molecular weight by Western-blot: 14 kDa

Preparation

This antibody was obtained by immunization of PCBD2 recombinant protein (Accession Number: NM_032151). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:113684(PCBD2 antibody) at dilution of 1:500