

ABCB6 antibody

Catalog Number: 107684

Product name

ABCB6 antibody

IHC: 1:20-1:200

IF: 1:20-1:200

Specificity

Human; other species not tested.

Antibody description

ABCB6 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive IF detected in MCF-7 cells. Positive WB detected in Raji cells. Observed molecular weight by Western-blot: 170 kDa

Preparation

This antibody was obtained by immunization of ABCB6 recombinant protein (Accession Number: NM_005689). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:107684(ABCB6 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:107684(ABCB6 antibody) at dilution of 1:50 (under 10x lens)



Immunofluorescent analysis of MCF-7 cells using

ABCB6 antibody



Catalog Number: 107684

Catalog No:107684(ABCB6 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)