

Human PRSS22/BSSP4/PRSS26/SP001LA (His tag) recombinant protein



Catalog Number: 519684

General Information

Gene Name Synonym

Brain-Specific Serine Protease 4, BSSP-4, Serine Protease 22, Serine Protease 26, Trypsin Epsilon, PRSS22, BSSP4, PRSS26

Protein Construction

A DNA sequence encoding the Human PRSS22/BSSP4/PRSS26/SP001LA Ala33-Ser317 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM HAc-NaAc, 150mM NaCl, 10% Glycerol, pH 4.5.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human PRSS22/BSSP4/PRSS26/SP001LA (His tag) recombinant protein