

ROMO1 antibody

Catalog Number: 114784

Product name

ROMO1 antibody

Specificity

Human; other species not tested.

Antibody description

ROMO1 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue, human gliomas tissue.

Preparation

This antibody was obtained by immunization of ROMO1 recombinant protein (Accession Number: XM_017027678). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC

Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations

Immunohistochemistry of paraffin-embedded human lung cancer slide using Catalog No:114784(ROMO1 Antibody) at dilution of 1:50