

# SYNJ1 antibody



Catalog Number: 115784

## Product name

SYNJ1 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

SYNJ1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, NIH/3T3 cells. Positive IF detected in NIH/3T3 cells. Observed molecular weight by Western-blot: 140 kDa

## Preparation

This antibody was obtained by immunization of SYNJ1 recombinant protein (Accession Number: BC136603). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:20-1:200

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115784(SYNJ1 Antibody) at dilution of 1:4000



Immunofluorescent analysis of (-20oc Ethanol) fixed NIH/3T3 cells using Catalog No:115784(SYNJ1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)