

Anti-LILRB1 Rabbit antibody



Catalog Number: 167784

Product name

Anti-LILRB1 Rabbit antibody

Specificity

Human

Antibody description

LILRB1 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human LILRB1

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/1000

ICC/IF: 1/20

FC: 1/50

IP: 1/20

Validations



Western blot detection of LILRB1 in HL-60,U251 cell lysates using LILRB1 Rabbit mAb(1:1000 diluted).Predicted band size:71kDa.Observed band size:110kDa.