

Cytochrome b5 antibody



Catalog Number: 109784

Product name

Cytochrome b5 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Cytochrome b5 Rabbit Polyclonal antibody.
Positive WB detected in human liver tissue.
Positive IHC detected in human pancreas cancer.
Observed molecular weight by Western-blot: 15kd

Preparation

This antibody was obtained by immunization of Cyt^ochrome b5 recombinant protein (Accession Number: NM_148923). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:109784(CYB5A antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:109784(CYB5A antibody) at dilution of 1:50 (under 10x lens)