

DNAJB3 antibody

Catalog Number: 110094

Product name

DNAJB3 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

DNAJB3 Rabbit Polyclonal antibody. Positive IHC detected in human testis. Positive WB detected in mouse testis tissue, Jurkat cells, PC-3 cells. Observed molecular weight by Western-blot: 28-30kd

Preparation

This antibody was obtained by immunization of DNAJB3 recombinant protein (Accession Number: NM_001001394). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200


Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:110094(DNAJB3 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human testis using Catalog No:110094(DNAJB3 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human testis using Catalog No:110094(DNAJB3 antibody) at dilution of 1:100 (under 40x lens)