

Rabbit IL12B/IL-12B (His Tag) recombinant protein



Catalog Number: 501094

General Information

Gene Name Synonym

Cytotoxic lymphocyte maturation factor 40 kDa subunit; IL-12 subunit p40

Protein Construction

A DNA sequence encoding the rabbit IL12B (XP_002710393.1) (Met1-Asn324) was expressed with a polyhistidine tag at the C-terminus.

Organism

Rabbit

Expression Host

Human Cells

QC Testing

Purity

> (29.1+65.8) as determined by SDS-PAGE.

Endotoxin

< 1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ile 23

Molecular Mass

The recombinant rabbit IL12B consists 313 amino acids and predicts a molecular mass of 36.1 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Rabbit IL12B / IL-12B Protein (His Tag)