

Nectin 2 antibody



Catalog Number: 113094

Product name

Nectin 2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Nectin 2 Rabbit Polyclonal antibody. Positive FC detected in MCF-7 cells. Positive WB detected in mouse kidney tissue, 3T3-L1 cells, mouse heart tissue. Observed molecular weight by Western-blot: 64 kDa

Preparation

This antibody was obtained by immunization of Nectin 2 recombinant protein (Accession Number: NM_002856). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, FC

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:113094(Nectin 2 antibody) at dilution of 1:300



1X10⁶ MCF-7 cells were stained with 0.2ug Nectin 2 antibody (Catalog No:113094, red) and control antibody (blue). Fixed with 4% PFA blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.