

# Anti-MRPL3 antibody



Catalog Number: 176094

## Product name

Polyclonal

Anti-MRPL3 antibody

## Ig Type

## Specificity

Rabbit IgG

Human

## Applications

## Antibody description

WB

Rabbit polyclonal antibody to MRPL3

## Dilutions

## Preparation


WB: 1:500-1:2000

This antigen of this antibody was synthesized peptide derived from mrp-l3 . at aa range: 20-100

## Validations

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.



## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Fig1.; Western Blot analysis of various cells using MRP-L3 Polyclonal Antibody

## Clonality