

TMEM120A antibody



Catalog Number: 116094

Product name

TMEM120A antibody

Immunogen

[Human TMEM120A Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TMEM120A Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 45 kDa

Preparation

This antibody was obtained by immunization of TMEM120A recombinant protein (Accession Number: BC029487). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:10-1:100

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116094(TMEM120A Antibody) at dilution of 1:300



Immunofluorescent analysis of HeLa cells using Catalog No:116094(TMEM120A Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)