

Phospho-STK11-T189 pAb



Catalog Number: 166494

Product name

Phospho-STK11-T189 pAb

Specificity

Human

Antibody description

Polyclonal antibody to Phospho-STK11-T189

Preparation

Antigen: A synthetic phosphorylated peptide around T189 of human STK11 (NP_000446.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - Phospho-STK11-T189 pAb

Western blot analysis of extracts of A431 cells, using Phospho-STK11-T189 antibody at 1:1000 dilution. A431 cells were treated by 10% FBS for after serum-starvation overnight. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% BSA.