

TULP4 antibody



Catalog Number: 116494

Product name

TULP4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TULP4 Rabbit Polyclonal antibody. Positive WB detected in human skeletal muscle tissue, human brain tissue. Positive IHC detected in human kidney tissue. Observed molecular weight by Western-blot: 169kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_020245). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:

WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:116494(TULP4 antibody) at dilution of 1:800