

# Anti-HBEGF/DTR Rabbit antibody



Catalog Number: 167494

## Product name

Anti-HBEGF/DTR Rabbit antibody

## Specificity

Human, Rat

## Antibody description

HBEGF/DTR Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human HBEGF

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC-P

## Dilutions

WB: 1/2000-1/10000

IHC: 1/20-1/100

## Validations



Western blot detection of HBEGF in K562, C6, 3T3, HeLa cell lysates using HBEGF Rabbit pAb(1:1000 diluted). Predicted band size: 23kDa. Observed band size: 23kDa.