

# Anti-Collagen I Rabbit antibody



Catalog Number: 166894

## Product name

Anti-Collagen I Rabbit antibody

## Specificity

Human, Cow

## Antibody description

Collagen I Rabbit polyclonal antibody

## Preparation

Antigen: A synthesized peptide derived from human Collagen I

## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC-P

## Dilutions


WB: 1/500-1/2000

IHC: 1/50-1/200

## Validations



Western blot analysis of Collagen I in Human stomach tissue lysate.



Immunohistochemical analysis of paraffin-embedded human kidney, using Collagen I Antibody.