

# Anti-Peroxiredoxin 3/PRDX3 Rabbit antibody



Catalog Number: 167894

## Product name

Rabbit IgG

Anti-Peroxiredoxin 3/PRDX3 Rabbit antibody

## Applications

## Specificity

WB, IHC-P, ICC/IF, FC, IP

Human

## Dilutions

## Antibody description

WB: 1/2000-1/10000

Peroxiredoxin 3/PRDX3 Rabbit monoclonal antibody

IHC: 1/50

ICC/IF: 1/50

## Preparation

FC: 1/20

Antigen: A synthetic peptide of human Peroxiredoxin 3/PRDX3

IP: 1/20-1/100

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Validations


## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type



Western blot detection of Peroxiredoxin 3/PRDX3 in HeLa,HL-60, using Peroxiredoxin 3/PRDX3 Rabbit mAb(1:1000 diluted)